

Title	Rabies Specimen Collection and Transport
Short Name	Rabies
Specimen Requirements	<ol style="list-style-type: none"> 1. Since brain tissue is examined for the diagnosis of rabies, only the animal head will be accepted for diagnostic purposes. Animals should be euthanized in a manner that will not damage the brain. The neck should be severed at the midpoint between the base of the skull and shoulders. 2. For bats, the entire dead animal should be submitted. 3. Only fresh, non-fixed tissue is acceptable. <p><u>The ISDH Rabies Laboratory requires that animal specimens meet the following criteria before accepting and testing specimens:</u></p> <ol style="list-style-type: none"> 1. Any animal that has bitten or otherwise exposed a human, pet or livestock: <ol style="list-style-type: none"> a. Observation of dogs, cats, and ferrets for 10 days is required unless the animal is unable to be observed due to injuries, illness at the time of the incident, or becomes ill during the observation period. b. Submission of a dog, cat, or ferret under other circumstances requires prior notification and approval of the ISDH Veterinary Epidemiologist. c. Any other domestic or wild animal that is suspected of being rabid, has bitten or otherwise exposed a human, pet, or livestock should be submitted as soon as possible. 2. Bats that have been involved in biting or non-bite exposures with humans or animals. <ol style="list-style-type: none"> a. Bats found in living space of human dwellings where people have been sleeping (an individual awakens to find a bat in the room, or an adult observes a bat in a room with an unattended child, or a mentally disabled or intoxicated person) may be submitted even though actual exposure cannot be confirmed. b. Bats found sick, dying, or dead outdoors where there is no known or suspected exposure to humans or domestic animals should not be submitted. 3. Terrestrial species that are known to be rabies vectors in Indiana and other states (skunks, foxes, raccoons, coyotes) whose actions suggest abnormal behavior or central nervous system (CNS) disease may be submitted whether a bite or non-bite exposure to a human or animal has occurred. These animals will be tested for surveillance purposes. 4. Squirrels, mice, hamsters, gerbils, rats, or other animals that have a very low or non-incidence of rabies will not be tested unless there is a human bite or other exposure involved and with prior approval of the ISDH Laboratory Director or Veterinary Epidemiologist. 5. Brain material from animals with CNS disease symptoms may be submitted by veterinarians or veterinary diagnostic laboratories as part of their diagnostic process, regardless of bite or other exposure status.
Sampling Materials	<ol style="list-style-type: none"> 1. Leak-proof containers (e.g., heavy-duty plastic bag). 2. Cold packs. 3. Zip-lock bag for form. 4. Insulated shipping boxes/containers.
Procedural Notes	<ol style="list-style-type: none"> 1. All specimens should be stored at refrigerator temperature (2-8°C) until shipped. Freezing of specimens is discouraged because it may delay testing since it may take up to 24 hours to thaw; it also damages the tissue. 2. If multiple specimens are shipped in the same container, be sure each specimen is labeled so it can be matched with the appropriate submission information. 3. Specimens should be submitted by LimsNet, an online system that will make results available as PDF files the minute they are released at the lab. 4. To get a free LimsNet account established at your facility for electronic submission and results reporting, call the help desk at (888) 535-0011 or email to LimsAppSupport@isdh.in.gov. 5. Alternatively, a totally completed Rabies Examination Form must accompany each specimen.
Shipping Instructions	<ol style="list-style-type: none"> 1. Place animal head specimens or whole bat specimens in a leak proof container (e.g., double bag using heavy-duty plastic bags or a metal can with lid) and seal tightly. 2. Place the leak proof container in an insulated shipping carton (e.g., styrofoam shipper) with sufficient cold packs to maintain refrigeration temperature during shipment. Do not use wet ice or packing peanuts.

3. Place the completed examination form in a zip-lock bag before placing it in the shipping container.
4. Place the insulated shipping carton into a cardboard shipping box, close lid, and seal.
5. Shipment via courier service is usually the most rapid method of delivery. The US Postal Service will not accept this type of specimen for mailing. At this time, United Parcel Service (UPS) will accept specimens for delivery to the ISDH Labs if properly packaged. Consult with your local UPS representative to obtain the most current shipping/packaging requirements for UPS.
 - a. Ship to:
Indiana State Department of Health Laboratories
Rabies/Virology Laboratory
550 W. 16th Street, Suite B
Indianapolis, IN 46202